



68 sec. video

Product Catalogue




oscident[©]

STRIPS

Strips made of medical steel, double-sided diamond coated, sterilisable, packed in a transparent single blister.

Article	Description	REF-No.	Price
● Strip (15µm)	Polisher, extra fine	1015	€20
● Strip (25µm)	Strip perforated, fine	1025	€20
● Strip (40µm)	Strip perforated, medium	1040	€20
● Strip (60µm)	Strip perforated, coarse	1060	€20
● Strip (90µm)	Strip perforated, extra coarse	1090	€20
● Strip (15µm)	Opener extra fine	1115	€22
● Strip (40µm)	Opener medium	1140	€22



STRIP SETS IN ALUMINIUM CASE

Individual selection of OPS and StarterSet, available in aluminium case or transparent blister.

Article (amount)	REF-No.	Price
OnePatientSet (3) in alu-case ● ● ●	1003	€69
StarterSetStrips (6) incl. measuring gauge in alu-case ● ● ● ● ● ●	1006	€129
StarterSetASR (6) incl. measuring gauge in alu-case		
like REF-No. 1006, incl. 1 Sirona T3LineEVA 11	1009A	€759
like REF-No. 1006, incl. 1 W&H WG-69 A	1009B	€759
like REF-No. 1006, incl. 1 NSK Ti-Max X55	1009C	€599



1003



1006

OSCILLATING EVA CONTRA ANGLE HANDPIECES

The oscillating EVA contra angle handpieces are easy to work with, minimally-invasive, and highly efficient.

Article	Description	REF-No.	Price
Sirona T1LineEVA 11L	2.4:1 Ratio, nozzle spray, stroke movement 1.1 mm, with fiber optics	1110	€ 789
Sirona T3LineEVA 11	2.4:1 Ratio, nozzle spray, stroke movement 1.1 mm, without fiber optics	1112	€ 659
NSK Ti-Max X55L	4:1 Ratio, nozzle spray, stroke movement 1.4 mm, with fiber optics	1210	€ 679
NSK Ti-Max X55	4:1 Ratio, nozzle spray, stroke movement 1.4 mm, without fiber optics	1211	€ 499
W&H WG-69 LT	2:1 Ratio green, nozzle spray, stroke movement 0.9 mm, with fiber optics	1252	€ 897
W&H WG-69 A	2:1 Ratio green, nozzle spray, stroke movement 0.9 mm, without fiber optics	1253	€ 623



PROSTHETIC SCREW DRIVER

cordless, time saving, versatile, easy to use, precisely adjustable torque.

Article	Description	REF-No.	Price
NSK ISD-900	Reliable and durable handpiece with excellent handling control, accoustic signal for reverse running, easy replaceable battery, LCD-operating-panel, easy to read and to use, torque control calibration system (TCS)	1320	€1.395
W&H IA-400	Precise torque control adjustable from 8 – 40 Ncm in 1 Ncm steps, with one finger-switch, independent shut-off while reaching the torque, time saving due to constant speed, wireless working, lifetime lubricated, reduction of ratchet systems	1321	€1.499



oscident[®]

oscident GmbH

Bommersheimerweg 11a
61348 Bad Homburg v.d.H.
Germany

Phone +49 (0)6172 2663956
Email info@oscident.de
Internet www.oscident.de

Misprints and errors excepted.
All prices exclude VAT and shipping.