

Strip band sizes

Article	Band size*	Grain size	Description	Double-sided REF-NR.	One-sided L REF-NR.	One-sided R REF-NR.
 Polisher	0.10mm	15µm	extra fine	1015D	1015L	1015R
 Strip	0.15mm	25µm	fine	1025D	1025L	1025R
 Strip	0.20mm	40µm	medium	1040D	1040L	1040R
 Strip	0.25mm	60µm	coarse	1060D	1060L	1060R
 Strip	0.30mm	90µm	extra coarse	1090D	1090L	1090R
 Strip	0.40mm	40µm	medium	1145D		
 Strip	0.50mm	90µm	extra coarse	1150D		
 Opener	0.10mm	15µm	fine	1115D		
 Opener	0.20mm	40µm	medium	1140D		



*approx. band size during first use of the instrument. Always check IPR space with a measuring gauge.
All strips available in packages of 1pcs and 5pcs. Prices are quoted in EUR and do not include VAT or shipping.